

The LumiR User Guide

Renan Sauteraud*

March 18, 2013

Contents

1	Introduction	1
1.1	Technology	1
1.2	LumiR	1
2	blum	1
3	bsum	10

1 Introduction

1.1 Technology

Luminex xMAP technology is a multiplex assay using flow cytometry to concurrently measure up to 500 analytes in a single reaction volume.

1.2 LumiR

LumiR is an R package that provides data structure and functions to read, store and analyze Luminex xMAP experiment.

As with any R package, it should first be loaded in the session

```
> library(LumiR)
```

2 blum

The first step is to read the data into R using the function `read.experiment`. Assuming the folder to read meets all the requirements, the function will automatically detect what software has been used for the acquisition. `read.experiment` returns an object of class `blum`. This object contains the bead level data and all the information extracted from the mapping files.

```
> XP.path <- system.file("extdata", package = "LumiR")
> bl <- read.experiment(XP.path)
```

*rsautera@fhcrc.org

phenotype mapping file missing.
Some information may be missing.

0%

0%

0%

89%

0%

0%

0%

0%

0%

0%

0%

0%

0%

68%

94%

0%

0%

0%

0%

98%

0%

0%

0%

0%

0%

```
> bl
```

An object of class blum with 27 analytes:

CHEX.1 CHEX.2 CHEX.3 ... IL.7 IL.8 IP10

885047 measures of expression in 176 wells, on 2 plates.

And slots: phenoData featureData exprs

LumiR makes use of the `ggplot2` package to quickly visualize the design of an experiment. The `plot_layout` method can plot phenotype data for the selected plate.

```
> plot_layout(bl, plate = "plate1", fill = "sample_type")
```

Here, we use colors to plot the `sample_type` for each well of `plate1`.

This method uses `ggplot2`'s `geom_polygon` function, therefore, the same arguments can be used here.

3 bsum

The second step is to summarize the data using `slummarize`. It returns an object of class `slum`. It is used for 'classic' standardization methods as it only contains the median fluorescence intensity (MFI) for each well and not the raw bead events. For this step, the package needs a user provided 'layout.csv' file with at least the location and expected concentration of the standards. The file provided in the `extdata` folder of this package is an example of such file.

```
> sl <- slummarize(bl)
```

Since `slum` objects retain the phenotype information from the mapping file, they may be used as first argument in `plot_layout`.

`slummarize` does the standard curve fitting using a 5-Parameters Logistic curve by default. LumiR defines a new method to be used in `ggplot2` for standard curve visualisation.

In the following example, we melt the object into a big data.frame in order to display the standards FI on the same plot using `ggplot2::geom_point`.

```
> library(ggplot2)
> msl <- melt(sl)
> msl.ss <- subset(msl, tolower(sample_type) == "standard")
> ggplot(msl.ss, aes(color = plate), alpha = 0.5) + scale_x_log10() +
+   scale_y_log10() + facet_wrap(~analyte) + geom_sc(sl) + geom_point(aes(x = concentration,
+   y = mfi))
```

Here again, any argument will be passed to the underlying `ggplot2`'s method `geom_line`.